

# Export Basics

After you have iterated on your recipe and generated a result that is to your satisfaction, you can export the transformed data.

## Steps:

1. In the left nav bar, click the Jobs icon.
2. In the Jobs page, click the job identifier to open the job in the Job Details page.
3. Click the Output Destinations tab.

## Export by:

- **Direct file download:** Click the file to download. From its righthand context menu, select **Download result**.

**NOTE:** Some file types cannot be downloaded.

- **Create new dataset:** You can create a new dataset from a generated output. Click the file. From its righthand context menu, select **Create imported dataset**.
- **Publish:** If Trifacta® has been integrated with an external datastore, you can publish your results to a designated target. Click **Publish**.